



Protocol for Testing Metal Cookware for Lead and Other Toxic Metals

Contents

- Introduction 2
- Acknowledgements 2
- Disclaimer 3
- Acronyms 3
- 1. Screening Cookware: Selection of Items 4
- 2. Preliminary Data 4
 - Characterizing the Cookware 4
- 3. XRF Preparation 7
 - 3.1. Health and Safety 7
 - 3.2. XRF Calibration 7
- 4. Sample Preparation 9
- 5. XRF Operation 10
- 6. XRF Data Documentation 10
- 7. Data Back-Up 13
- Data Documentation 14
- Appendix A: XRF Health and Safety 15

Introduction

The main objective of this protocol is to present a step-by-step process for screening metal cookware for lead and other toxic metals (e.g. cadmium, arsenic, etc.) using an XRF device. The data entry procedures described in this document assume direct entry of information into a [series of spreadsheets](#). However, [templates are also available](#) to handwrite results on paper. The handwritten information should then be entered into the spreadsheets.

Acknowledgements

The following people reviewed and provided comments on the document, and their contributions are gratefully acknowledged:

Gordon Binkhorst | Pure Earth (October, 2024)
Katie Fellows | King County (April, 2024)
Richard Fuller | Pure Earth (May, 2024)
Angela Mathee | SAMRC (April, 2024)
Kate Porterfield | Pure Earth (October, 2024)
Mahbubur Rahman | iccdr,b (May, 2024)
Talia Reiss | Pure Earth (July 2024)
Alfonso Rodriguez | Pure Earth (October, 2024)
Renee Street | SAMRC (April, 2024)
Jesmin Sultana | iccdr,b (May, 2024)
Budi Susilorini | Pure Earth (May, 2024)
Stephen G. Whittaker, Ph.D. | Pure Earth (October, 2024)

Disclaimer

This protocol is designed for use in specific projects and may not be universally applicable. It should be adapted or modified only with the guidance of relevant experts to ensure it meets the unique needs of each project. The creators of this protocol assume no responsibility for its misuse or for any outcomes resulting from its application beyond its intended scope.

Acronyms

HBA	Home-based assessment
RMS	Rapid Market Screening
XRF	X-ray Fluorescence

1. Screening Cookware: Selection of Items

When analyzing for a Rapid Market Screening (RMS) or Home-Based Assessment (HBA), please reference [RMS](#) and [HBA](#) protocols regarding the criteria for measuring the total lead content of metal cookware. The main points are summarized below:

Select the most common items used for cooking (e.g., pertinent to a home, school, marketplace) in consideration of the project objectives (e.g., aluminum cookware versus brass or stainless steel).

2. Preliminary Data

Characterizing the Cookware

This section is designed to record the general characteristics of the cookware to be tested. Please complete as many of the following data fields as possible in the [provided template](#):

Table 1: Cookware Characteristics

Variable	Data Entry	Definition
Investigator	Free-form	Full name of investigator
OrganizationTesting	Free-form	Organization conducting the analysis
CountryAcquiredFromCode	Dropdown	Alpha-3 code for the country where the cookware was acquired from
CityAcquiredFromCode	Free-form	3-digit abbreviation for the city where the cookware was acquired from (UNK if unknown)

Variable	Data Entry	Definition
CookpotNumber	Dropdown	Assigned cookpot number (001, 002, 003, etc.)
Cookpot_ID	Autogenerated	<u>Autogenerated</u> Cookpot ID from the combination of the previous three fields (for linking to XRF and leachate data)
ManufType	Dropdown	Was the cookpot manufactured by casting molten metal or forged (i.e., shaping of ingots)? UNK if unknown.
CookpotShape	Free-form	Describe the shape of the cookpot (i.e., flat or round bottom, spherical or typical frying pan shape, wok, etc.)
CookpotWidth	Free-form	Average internal width of cookpot (in centimeters)
CookpotDepth	Free-form	Average internal depth of cookpot (in centimeters)
InsideDiffOutside	Dropdown	Is the inside of the pot visibly different from the outside? Yes or No
Photos	Free-form	Photo filenames/links (for pots, components, and labels). For example: 2345.jpg
Manufacturer	Free-form	Manufacturer of cookpot. UNK if unknown.
MakeModel	Free-form	Make/Model of cookpot. UNK if unknown.

Variable	Data Entry	Definition
Volume	Free-form	Volume of pot in quarts or liters (specify). UNK if unknown.
CookpotType	Dropdown	Type of cookpot: pressure cooker, fry pan, sauté pan, saucepan, wok, brazier, stock pot, fryer pot, steamer, caldero, double boiler, kadai, uruli, tadka pan, appam pan, idli maker, tope, handi, dadesen, other
OtherDesc	Free-form	"Other" cookpot description
Description	Free-form	Description of item from the online listing, packaging, or a description of donated items.
Condition	Dropdown	New (purchased locally or online) or Used (donated by the community) or Unknown/UNK
Metal	Dropdown	Primary cookpot metal: Aluminum, Stainless Steel, Hindalium (if noted by the seller of the item), Brass, Other, or Unknown (UNK)
Coating	Dropdown	Yes for an observed non-stick coating, No for no observed coating, UNK for Unknown
Anodized	Dropdown	Yes for manufacturer-noted anodization, No for no mention of anodization, UNK for Unknown
Certification	Free-Form	Certification listed on packaging or item (e.g., NSF or UL).

Variable	Data Entry	Definition
Notes	Free-Form	Additional narrative
Definitions: <i>Autogenerated</i> fields are populated automatically based on data entered in other fields <i>Dropdown</i> fields restrict data entry to defined choices from a dropdown menu <i>Free-form</i> means that the analyst can enter unstructured text, without specific formatting		

3. XRF Preparation

3.1. Health and Safety

Before operating an XRF analyzer, read and follow the Health and Safety section in the Detailed XRF Guidance in [Appendix A](#) of this section.

3.2. XRF Calibration

- Inspect the measurement window for damage (rips) or dirt/debris – circled in red below. Replace if needed.



Figure 1. XRF Measurement Window

- Conduct a **calibration verification**. Calibration procedure is determined by the XRF analyzer model in use:
- **If using a NITON XRF Machine:**
 - Press the “System Check” icon to initiate the system check process.
 - Once the instrument completes the check, press the “Sample Type” icon, followed by the “Metals” icon, and then the “General Metals” icon.
 - Measure the alloy standard sample: Test the factory supplied alloy standard (or other approved standard).
 - If the standard is correctly identified, and all major elements read within calculated acceptance limits (see factory QC readings), testing may commence.
 - If the standard is not within acceptance limits, turn the device off and on again, and then conduct a second systems check.
 - After the instrument completes the system check, return to “Test” and re-measure the alloy standard as directed above. If it reads correctly, testing may proceed.
 - **Significant deviations from the expected values may indicate a need for recalibration or adjustment of the instrument. Contact your Thermo Scientific Niton Service center for assistance and do NOT proceed with testing.**
- **If using an OLYMPUS XRF Machine:**
 - Calibration checks (Cal Checks) can be performed in two ways:
 - The docking station provides an automatic Cal Check when inserted.
 - When in the field or away from the docking station, utilize the standardization coupon included in the XRF kit to perform a Cal Check.
 - To perform a manual Cal Check:
 - Navigate to the Test Setup screen
 - Place the supplied Cal Check alloy standard (316 stainless steel) on a flat surface
 - Ensure the analyzer measurement window is flush over the coupon
 - Tap the Cal Check button
 - The message “Cal Check - Passed” indicates that the analyzer is ready for sample testing
 - If the Cal Check fails:
 - Ensure that the alloy standard is positioned correctly beneath the measurement window

- Confirm that the X-ray indicator blinks during the procedure
 - Restart the analyzer
 - Retry the Cal Check procedure
 - **If the Cal Check fails repeatedly, contact Olympus customer service or a local distributor**
- If you are not using a check sample that was included in the XRF kits, refer to the following standards for each sample type or matrix:

Table 2: XRF Standards by Metal Sample Type

Matrix	Expected Pb conc (ppm)	Lower acceptable range (~25%)	Upper acceptable range (+25%)
MBH-1611X SAC305Q (“Tin Check”)	1000	750	1250
MBH Check Sample R 180-696 (Alternate Tin Standard)	1200	900	1500

- The frequency of calibration checks depends on several factors, including the instrument usage, environmental conditions, and regulatory requirements. As a general guideline for analyzing metal cookware, perform the calibration check on the XRF machine as follows:
 - Daily or before each analysis session;
 - At the beginning session;
 - At the end of the session;
 - When the instrument has been switched off for 30 minutes;
 - After analyzing every fifth piece of cookware.
- Record lead measurement result for each calibration standard used.
- Do not proceed if pertinent calibration standard checks are outside of the acceptance limits of the device.

4. Sample Preparation

The following preparation guidelines should be followed for metal cookware:

- When accessible, samples should be taken from the interior of the cookware, where food would be in contact.
- For curved or irregularly shaped items, analyze the flattest surfaces.
- Retain all samples (properly labeled and secured) for possible additional screening or laboratory testing.

5. XRF Operation

- Samples should be placed on a hard surface, and **not held in your hand**, when analyzing.
- Select TEST ALL mode on the XRF device. (if applicable)
- Take an XRF reading of the table or surface you will be using to analyze the samples to ensure it does not contain lead, as this could interfere with the readings. If lead is detected, find a new surface.
- Take the following XRF readings; NOTE: XRF readings should be about 30 to 60 seconds long depending on the model used:
 - 6 readings of the interior sides of the cookware, where food would be in contact
 - 6 readings of the interior bottom of the cookware
 - 6 readings of the interior cover/lid of the cookware (if applicable)
 - If applicable, the following readings should be included:
 - 2 reading of each rivet, knob, or other external features (e.g., valves for pressure cookers).
 - Record what part of the item each reading corresponds to.
- If the interior of the cookware is not accessible (e.g., too small), but it is apparent that the interior and exterior materials are similar in makeup, take XRF readings as follows:
 - 6 readings of the exterior sides of the cookware
 - 6 readings of the exterior bottom of the cookware
 - 6 readings of the exterior top of the cookware (if applicable)
 - If applicable, the following readings should be included:
 - 2 reading of each rivet, knob, or other external features of the cookware

6. XRF Data Documentation

In addition to recording the standard cookware information (see Table 1), please complete as many of the following data fields from Table 3 as possible in the [provided template](#).

Table 3: XRF Data Fields

Variable	Data entry	Definition
User	Free-Form	Full name of XRF user
OrganizationTesting	Free-form	Organization conducting the analysis
CountryTested	Dropdown	Country where the cookware was tested: Alpha-3 code
CityTested	Free-Form	3-digit abbreviation for the city where the cookware was tested
XRF_Model	Free-Form	Make and model of XRF analyzer
XRF_Mode	Free-Form	XRF settings used
Date	Format enforced	Date of XRF measurement
Time	Format enforced	Time of XRF measurement.
Cookpot_ID	Free-form	From Cookpot ID generated in General Dictionary
TestedComponent Designation	Dropdown	Letter assigned to the cookpot component being tested (a,b,c,d,e, etc. If >26 components, start with aa, then ab, ac, etc.)

Variable	Data entry	Definition
SampleID	Autogenerated	Autogenerated Sample ID from the combination of TestedComponentDesignation and Cookpot_ID
XRF_Run_No	Free-Form	XRF run number
InsideXRFMeasure	Dropdown	Is the inside of the pot (i.e., the cooking surface) accessible by the XRF? Yes or No
ReadingLocation	Dropdown	outside lid; inside lid; outside pot base; inside pot base; outside pot wall; inside pot wall; handle; inside steamer insert; outside steamer insert; rivet; vent pipe; other
ComponentMaterialTested	Dropdown	Metal, plastic, glass, or rubber, or other, as determined upon inspection.
FoodContact	Dropdown	Yes/No whether the component could contribute to lead in food or leachate.
Pb	Free-Form	Concentration detected by the XRF. If non-detected, enter the "<" symbol and enter the error term in the field "Error". Do not enter "0"
Units	Dropdown	XRF concentration units (ppm or %)
Error	Free-form	Error for the measurement, as reported by the XRF
ErrorTerm	Dropdown	Whether the instrument provides the error term as two or three standard deviations.

Variable	Data entry	Definition
XRF_Data_File	Free-form	File containing the XRF reading on the instrument
Notes	Free-form	Additional narrative
<p>Definitions:</p> <p><i>Autogenerated</i> fields are populated automatically based on data entered in other fields</p> <p><i>Dropdown</i> fields restrict data entry to defined choices from a dropdown menu</p> <p><i>Format enforced</i> fields restrict data entry to specific formats</p> <p><i>Free-form</i> means that the analyst can enter unstructured text, without specific formatting</p>		

7. Data Back-Up

- At the end of each day when XRF readings are taken, download and retain raw data files as .csv or .xls.
- NOTE: Make sure that each XRF reading in the data file corresponds to each sample reading through the XRF Reading Log; otherwise, the data file is not useful.
- In the file name, include the date that the samples were analyzed as DD-MM-YYYY.

If the average or median of the internal surface XRF readings exceeds the reference value (100 ppm), a leachate test of this item may be warranted depending on the objectives of the particular study. For this follow the specific protocol [here](#).

Data Documentation

A downloadable Excel version of the sample data collection table can be found [here](#).

Appendix A: XRF Health and Safety

Correct operation of the XRF is critical for keeping investigators safe and for collecting accurate data.

All investigators operating the XRF should be familiar with the safety information presented in *Chapter 2: Using your analyzer of Niton XL3 Analyzer User's Guide*. The use of XRF is restricted to those people who have taken one of the training courses from one of our providers

XRF (X-ray fluorescence) is a non-destructive analytical technique used to determine the elemental composition of materials. XRF analyzers determine the chemistry of a sample by measuring the fluorescent (or secondary) X-ray emitted from a sample when it is excited by a primary X-ray source. Each of the elements present in a sample produces a set of characteristic fluorescent X-rays (“a fingerprint”) that is unique for that specific element. If you are interested in reading more about how the XRF works, read ThermoFisher’s [XRF Technology in the Field](#).

Primary radiation is radiation that is produced by the analyzer and emitted out through the measurement window.

- **Always treat radiation with respect.**
- **Do not hold your analyzer near the measurement window during testing. Never point your analyzer at yourself or anyone else when the shutter is open**
- There should always be a sample in contact with the measurement window when the x-ray tube is on.
- The sample will absorb most of the primary-beam radiation unless it is smaller than the instrument's measurement window or of low density and/or thickness. Caution should be taken when analyzing samples that are small, thin, and/or low in density as they may allow much more of the primary beam to escape.

The primary beam is a directed beam out of the front of the analyzer that can have high dose rates. The secondary beam, or scattered beam, has much lower dose rates.




Figure 2. Primary beam




Figure 3. Secondary (scattered) beam




Figure 4. X-ray beam indicator lights